


Human and Mouse MGC Collaboration Clones

Distributed through K.K. DNAFORM

Vector Map: pDONR-201

Vector Name: pDONR-201
 Source: Invitrogen
 Plasmid type: Gateway Donor Vector
 Selection marker: Kanamycin resistant; 50 µg/ml

Map



3661	tgtccttcta	gtgtagccgt	agttaggcca	ccacttcaag	aactctgtag	caccgcctac	3720
3721	atacctcgct	ctgctaatec	tgttaccagt	ggctgctgcc	agtggcgata	agtcgtgtct	3780
3781	taccgggttg	gactcaagac	gatagttacc	ggataaggcg	cagcggtcgg	gctgaacggg	3840
3841	gggttcgtgc	acacagccca	gcttggagcg	aacgacctac	accgaactga	gatacctaca	3900
3901	gcgtgagcta	tgagaaagcg	ccacgcttcc	cgaagggaga	aaggcggaca	ggtatccggt	3960
3961	aagcggcagg	gtcggaacag	gagagcgcac	gagggagctt	ccagggggaa	acgcctggta	4020
4021	tctttatagt	cctgtcgggt	ttcgccacct	ctgacttgag	cgtcgatttt	tgtgatgctc	4080
4081	gtcagggggg	cggagcctat	ggaaaaacgc	cagcaacgcg	gcctttttac	ggttcctggc	4140
4141	cttttgctgg	ccttttgctc	acatgttctt	tctgctgcta	tcccctgatt	ctgtggataa	4200
4201	ccgtattacc	gctagccagg	aagagtttgt	agaaacgcaa	aaaggccatc	cgtcaggatg	4260
4261	gccttctgct	tagtttgatg	cctggcagtt	tatggcgggc	gtcctgcccg	ccaccctccg	4320
4321	ggccgttgct	tcacaacggt	caaateccgct	cccggcggat	ttgtcctact	caggagagcg	4380
4381	ttcaccgaca	aacaacagat	aaaacgaaag	gcccgagtctt	ccgactgagc	ctttcgtttt	4440
4441	atattgatgcc	tggcagttcc	ctactctcgc	4470			

Contact Address

K.K. DNAFORM – *Precision Gene Technologies*
 Leading Venture Plaza-2
 75-1, Ono-cho, Tsurumi-ku, Yokohama City
 Kanagawa 230-0046
 Japan

Tel: +81-(0)45-510-0607
 Fax: +81-(0)45-510-0608
 e-mail: techinfo@dnaform.jp
 Homepage: <http://www.dnaform.jp>